

---

## Rat IgG1 Isotype Control SAFIRE Purified

Catalog Number :42813-25

RUO: For Research Use Only. Not for use in diagnostic procedures.

---

### Product Information

**Clone:** HRPN

**Format/Conjugate:** SAFIRE Purified

**Concentration:** 2.0 mg/mL

**Reactivity:** Mouse

**Laser:** Not Applicable

**Peak Emission:** Not Applicable

**Peak Excitation:** Not Applicable

**Filter:** Not Applicable

**Brightness (1=dim,5=brightest):** Not Applicable

**Isotype:** Rat IgG1

**Formulation:** Phosphate-buffered aqueous solution, pH7.2.

**Storage:** Product should be kept at 2-8°C.

**Applications:** FC

### Description

The HRPN monoclonal antibody is used as an Isotype Control immunoglobulin for rat IgG1 antibodies.

### Preparation & Storage

The product should be stored undiluted at 4°C. Do not freeze. The monoclonal antibody was purified utilizing affinity chromatography.

### Application Notes

The antibody has been analyzed for quality through the flow cytometric analysis of the relevant cell type. It is recommended that the isotype control be used at the same concentration as the primary antibody being tested.

### References