

Rat IgG2b Isotype Control SAFIRE Purified

Catalog Number :41823-25

RUO: For Research Use Only. Not for use in diagnostic procedures.

Product Information

Clone: LTF-2

Format/Conjugate: SAFIRE Purified

Concentration: 2.0 mg/mL

Reactivity: Mouse

Laser: Not Applicable

Peak Emission: Not Applicable

Peak Excitation: Not Applicable

Filter: Not Applicable

Brightness (1=dim,5=brightest): Not Applicable

Isotype: Rat IgG2b

Formulation: Phosphate-buffered aqueous solution, pH7.2.

Storage: Product should be kept at 2-8°C.

Applications: FC

Description

The LTF-2 monoclonal antibody is used as an Isotype Control immunoglobulin for rat IgG2b antibodies.

Preparation & Storage

The product should be stored undiluted at 4°C. Do not freeze. The monoclonal antibody was purified utilizing affinitychromatography.

Application Notes

The antibody has been analyzed for quality through the flow cytometric analysis of the relevant cell type. It is recommended that the isotype control be used at the same concentration as the primary antibody being tested.

References